

Introducing
Secur **MED**TM

MEDICATION MANAGEMENT CARTS

|O|R|G|A|N|I|Z|E|D|
for **SAFETY**

Model # 27034C-EL

Model # 27013C-PREL

Model # 27095P-EL

Model # 27090P-PREL

Model # 27085M-EL

Model # 27054M-MSEL

L LAKESIDE[®]
Ingenuity On Wheels[™]

1-800-558-8565 • www.eLakeside.com

Patient safety and efficient organization are essential requirements in the fast-paced world of healthcare. Your pharmacy staff needs to locate, assemble, and organize scheduled prescriptions on a constant basis. Nursing teams must then be able to transport and dispense the medications to patients with the highest level of security and efficiency. Lakeside Manufacturing understands these requirements and has developed its new line of SecurMED™ Medication Management Carts specifically to meet your needs.

SecurMED™

MEDICATION MANAGEMENT CARTS

|O|R|G|A|N|I|Z|E|D|

for
SAFETY

Model # 27013C-PREL



Model # 27034C-EL

SecurMED & JCAHO

As of January 1, 2004 the Joint Commission on Accreditation of Healthcare Organizations (JCAHO) adopted the following Medication Management Standards to improve patient safety. Lakeside Manufacturing's new SecurMED Medication Management Carts support these procedures.

JCAHO STANDARD MM.4.20:

Medications are Prepared Safely

JCAHO requires that all staff members who prepare medications use techniques to ensure accuracy and avoid contamination. All equipment used must be sterile and the work environment is to be organized, clean, and uncluttered, maintaining functionally separate areas for medication preparation.

SecurMED ADVANTAGES:

- > Large work surface top is easy to organize
- > Smooth surfaces are easy to clean and sanitize
- > Patient specific drawers are separated from work surface

JCAHO STANDARD MM.4.30

Medications are Appropriately Labeled

Medications must be labeled in a standardized manner according to organizational policy, applicable laws, current regulations, and standards of practice to minimize errors. Minimal standards include drug name, strength, amount, expiration date/time, and dispensing directions.

SecurMED™ ADVANTAGES:

- > Flexible label options for Med Bins and Cassettes
- > Patient specific trays and dividers for separating medications

JCAHO STANDARD MM.4.40:

Medications are Dispensed Safely

Medications must be dispensed in quantities to meet patient needs, but minimize diversion. Dispensing must occur in a timely manner to meet patient needs. Medications must be dispensed in the most ready-to-administer forms available from the manufacturer (i.e. unit dose formats).

SecurMED ADVANTAGES:

- > Unit dose dispensing in Med Bin Carts
- > Cart drawers and cassettes can be organized by patient, room, or floor

JCAHO STANDARD MM.4.50:

Medications Can Be Dispensed to Patients When the Pharmacy is Closed

The organization must have the ability to provide medications to meet urgent or emergency needs when the pharmacy is closed. JCAHO requires to process for providing after-hours medications be evaluated on an ongoing basis.

SecurMED ADVANTAGES:

- > Four locking options improve security and discourage tampering – Manual, Keypad, Proximity Reader Head and Mag Stripe Reader Head
- > Locking narcotics box keeps sensitive drugs secure
- > Custom labeling instructs night time nursing staff for patient medications

Cassette Dispensing

Cassette dispensing offers exceptional flexibility for patients in maximum bed environments.

- Medications can be hand-loaded or used with automated systems
- Cassettes fit in small or large platform carts
- Steel construction for durability and long-term use
- Cassette drawers have custom labels
- Cassettes slide out for easy loading in the pharmacy



Unit Dose Dispensing

Med Bin dispensing supports unit dose medications and creates visual medication separation.

- Drawers and trays designed to support unit dose medications
- Removable trays for easy loading
- Steel construction for durability and long-term use
- Ball bearing drawer slides provide smooth, efficient access



Punch Card Dispensing

Punch card dispensing allows maximum flexibility for 30 days of medication and patient supplies.

- Deep and wide drawers support 30-day medication punch cards
- Plenty of room for maintenance medications and supplies
- Wide handle drawer pulls for improved ergonomics and easy access
- Steel construction for durability and long-term use
- Ball bearing drawer slides provide smooth, efficient access



Compact Capacity

(Dimensions - 22½" x 23" x 49¼")



Security



Organization



Ergonomics

- > Space-efficient footprint fits into tight locations
- > Cassette option for multiple patients
- > Drawer options for patient medications and supplies
- > Drawer depth - 6"
- > Smooth, easy to clean work surface top
- > Non-marring composite bumpers to protect walls and doorways
- > Steel construction for lifetime durability
- > Quiet, non-marking dual-wheel casters (2 locking; 2 directional)
- > Clean, contemporary healthcare appearance

|O|R|G|A|N|I|Z|E|D|
for
SAFETY



Model # 27034C-EL



Model # 27095P-EL



Intermediate Capacity

(Dimensions - 30 1/2" x 23" x 49 1/4")



Security



Organization



Ergonomics



Model # 27085M-EL

- > Unit dose configured for multiple patients
- > Separate plastic bins or drawer dividers per patient
- > Drawer options for patient medications and supplies
- > Drawer depth - 9"
- > Smooth, easy to clean work surface top
- > Non-marring composite bumpers to protect walls and doorways
- > Steel construction for lifetime durability
- > Quiet, non-marking dual-wheel casters (2 locking; 2 directional)
- > Clean, contemporary healthcare appearance

|O|R|G|A|N|I|Z|E|D|

for SAFETY



Full-Size Capacity

(Dimensions - 43½" x 23" x 49¼")



Security



Organization



Ergonomics

- > Deep and wide drawers support 30-day medication punch cards
- > Wide handle for improved ergonomics and easy access
- > Ball bearing drawer slides create smooth, efficient access
- > Separate plastic bins or drawer dividers per patient
- > Drawer depths – 3", 6" and 9"
- > Smooth, easy to clean work surface top
- > Non-marring composite bumpers to protect walls and doorways
- > Steel construction for lifetime durability
- > Quiet, non-marking dual-wheel casters (2 locking; 2 directional)
- > Clean, contemporary healthcare appearance



Model # 27013C-PREL



Model # 27054M-MSEL



Model # 27090P-PREL

|O|R|G|A|N|I|Z|E|D|
 for
SAFETY





CompX eLock's LockView® Software

LockView makes it easy to monitor audit trails for the last 1500 access attempts, manage users and eLocks, and assign time-based restrictions and dual credential access. Downloading and uploading information is easy. Simply attach one end of the RJ-11 cable to the back of the CompX eLock and the other end to a laptop computer – *there's no need for expensive hardwired systems!* LockView software allows for easy programming, audit trail retrieval and database maintenance. LockView is compatible with Windows™ 95 or newer.

- > Add, edit, view and delete users / supervisors and their credentials
- > Add, edit and delete CompX eLocks
- > Assign access rights between users / supervisors and CompX eLocks in the database
- > Download, view, save, print, archive and delete audit trails
- > Dual credential access available
- > Programmable lockout feature
- > Set real time clock
- > Time-based access restrictions
- > Create user groups
- > USB connection required

ORGANIZED
for
SAFETY

Model # 27054MS-EL



Security

Four Locking Options

- Key lock
- Electronic number pad
- Electronic card proximity
- Electronic card slide

Patient Labeling

- Computer-generated or hand written labels
- Plastic windows reduce label tampering
- Labels are easy to change

Narcotics Box

- Separate steel locking box for prescription narcotics
- Separate key for each box

Steel Construction

- Tamper-resistant metal construction with durable epoxy powder coat finish

Organization

Patient Specific Trays

- Molded plastic trays to organize patient medications
- Trays are removable for loading and cleaning

Drawer Dividers

- Dividers to organize patient medications
- Dividers can be arranged to different widths

Custom Labels

- Labels can be computer generated or hand written
- Labels are easy to remove and change
- Clear plastic window protects labels

Medication Dispensing

- Cassettes for high patient volume and automated systems
- Unit dose dividers and trays for maximum security
- Large open drawers for 30-day punch cards and supplies

Appearance and Cleaning

- Seamless work surface and smooth easy-to-clean finish
- Contemporary appearance fits healthcare environments

Ergonomics

Top Work Surface

- Large top for writing, organizing, and laptop use
- Three sides to keep supplies contained
- Smooth and easy to sanitize

Handle

- Plastic handle on top enables easy pushing

Casters

- High quality casters roll easily on hard floor or carpet
- Quiet and non-marking
- Wheel locks for stability in use

Composite Bumpers

- Bumpers on base protect walls and doorways
- Bumpers protect staff ankles from injury

Drawer Handles

- Finger slots allow for easy opening of drawers
- Full drawer width for easy access

SecurMED™ Medication Management Carts are available in three sizes with three dispensing options to fit the majority of acute and long term care operations.

Need a custom configuration? No problem! Lakeside Manufacturing's skilled design team can configure a cart to fit your specific application. The modular SecurMED platform allows you to mix and match components and options to fit your needs and budget. Best of all we offer 15-day turnaround on all custom orders!

Cassette Dispensing

MODEL	CART SIZE	DIMENSIONS			BIN WIDTH	# OF BINS	DRAWERS	# OF CASSETTES
		W	L	H				
27033C	Compact	22 ½"	23"	49 ¼"	3"	40	1 - 6" D	4
27034C	Compact	22 ½"	23"	49 ¼"	4"	32	1 - 6" D	4
27045C	Compact	22 ½"	23"	49 ¼"	5"	24	1 - 6" D	4
27048C	Compact	22 ½"	23"	49 ¼"	8"	16	1 - 6" D	4
27013C	Full-Size	43 ½"	23"	49 ¼"	3"	50	2 - 3" D; 1 - 6" D 2 - 9" D	5
27014C	Full-Size	43 ½"	23"	49 ¼"	4"	40	2 - 3" D; 1 - 6" D 2 - 9" D	5
27025C	Full-Size	43 ½"	23"	49 ¼"	5"	30	2 - 3" D; 1 - 6" D 2 - 9" D	5
27028C	Full-Size	43 ½"	23"	49 ¼"	8"	20	2 - 3" D; 1 - 6" D 2 - 9" D	5

Unit Dose Dispensing

MODEL	CART SIZE	DIMENSIONS			BIN WIDTH	# OF BINS	DRAWERS
		W	L	H			
27073M	Intermediate	30 ½"	23"	49 ¼"	3"	56	1 - 9" D 7 - 3" D
27074M	Intermediate	30 ½"	23"	49 ¼"	4"	42	1 - 9" D 7 - 3" D
27085M	Intermediate	30 ½"	23"	49 ¼"	5"	28	1 - 9" D 7 - 3" D
27088M	Intermediate	30 ½"	23"	49 ¼"	8"	21	1 - 9" D 7 - 3" D
27053M	Full-Size	43 ½"	23"	49 ¼"	3"	70	1 - 9" D 7 - 3" D
27054M	Full-Size	43 ½"	23"	49 ¼"	4"	56	1 - 9" D 7 - 3" D
27065M	Full-Size	43 ½"	23"	49 ¼"	5"	42	1 - 9" D 7 - 3" D
27068M	Full-Size	43 ½"	23"	49 ¼"	8"	32	1 - 9" D 7 - 3" D

Punch Card Dispensing

MODEL	CART SIZE	DIMENSIONS			# OF CARD HOLDER BRACKETS	DRAWERS
		W	L	H		
27095P	Intermediate	30 ½"	23"	49 ¼"	8	1 - 3" D 3 - 9" D
27090P	Full-Size	43 ½"	23"	49 ¼"	11	1 - 3" D 3 - 9" D

CHOOSE YOUR LOCKING SYSTEM

Add the following suffix to the above model numbers:

-M - Manual Key

-EL - Numeric keypad with override key

-PREL - Proximity with override key

-MSEL - Magnetic stripe with override key



LAKESIDE MANUFACTURING, INC.

4900 West Electric Ave.
West Milwaukee, WI 53219
U.S.A.

Toll Free: (800) 558-8565
Phone: (414) 902-6400
Fax: (414) 902-6446
e-mail: info@eLakeside.com
web: www.eLakeside.com

Call Lakeside Manufacturing at 1-800-558-8565 to start designing your SecurMED Medication Management Cart today!